



SERVER HD

MAX. LOAD CAPACITY 200 KG



SERVER HD

MAX. LOAD CAPACITY 200 KG



THE XL-BESTSELLER

SERVER HD was specially designed for powerful silhouettes. The rollator combines puristic design and clever ergonomic functions providing safety and comfort to users when driving, carrying or transporting the rollator. SERVER HD is available with popular TPE wheels or shock-absorbing soft wheels.

BENEFITS

- Large shopping bag for carrying up to 5 kg
- Ergonomic handles with locking function
- Handle height extremely convenient and easy to adapt
- Patented, red, "Click & Safe" safety catch
- Ingenious carrying handle: ergonomic handle grips under the mounting bars for lifting and carrying the rollator
- All-round safety reflectors available for shopping bag, frame and handles
- Large selection of accessories such as backrest straps, bags, tray, one-handed brake operation, slowdown brake and rollator umbrella

SIZES & COLOURS

Large: Anthrazit

Medium: Anthrazit



TECHNICAL DETAILS

Weight:	Max. load capacity	Seat height	Seat width	Width folded	Total width	Distance between handles	Handle heights from medium to large
Large 7.9 kg	bis 200 kg	Large 62 cm	55 cm	23 cm	69 cm	56 cm	66-102 cm
Medium 7.6 kg excl. shopping bag 340 grams		Medium 55 cm					
ZUBEHÖR: Rückengurt Klassik 8 cm breit					ZUBEHÖR: Rückengurt Komfort 20 cm breit		
86 cm Länge: Art. Nr. ACR01009XL					86 cm Länge: Art. Nr. ACR01009CX		
91 cm Länge: Art. Nr. ACR00910					91 cm Länge: Art. Nr. ACR00910C		

SPECIAL FEATURES



Sie können zwischen der TPE- und PUR-Softbereifung wählen




Inklusive Stockhalter für Ihre Gehhilfe




Ergo-Griffe mit Vertiefung als Unterarmablage



Bremshebel ohne viel Kraftaufwand bedienbar



Ergonomische Greifzone zum Heben und Tragen



Dank des „Click & Safe“ Verschlusses schließt der Rollator sicher

AUSWAHL AN ZUBEHÖR



Rückengurt



LED-Lampen



Transporttasche

Produktvideos auf



Ihr Rehasense Fachhändler