



EXPLORER

THE OUTDOOR ROLLATOR



EXPLORER

THE OUTDOOR ROLLATOR



MADE FOR OUTDOOR TERRAIN

The EXPLORER is the SUV among the rollators: this sporty model with many clever details was designed and manufactured especially for outdoor use. The flexi frame provides additional flexibility in rough terrain and makes the EXPLORER the perfect companion for your next outdoor adventure.

"My mother attends sport rehabilitation classes and loves spending time out in the fresh air so when I discovered the Explorer rollator I knew it would be perfect. She loves it and the independence she now has".

Nils from Bremen

FEATURES AT A GLANCE

- Large shopping bag for carrying up to 5 kg
- Ergonomic handles with locking function
- User friendly and easily adjustable handle height
- Patented, red, "Click & Safe" safety catch
- Ingenious carrying handle: ergonomic handle grips under the mounting bars for lifting and carrying the rollator
- All-round safety reflectors available for shopping bag, front and side of frame and handles

SIZES & COLOURS

Large, Medium: Brown



TECHNICAL DETAILS

Weight	Max. load capacity	Seat height	Seat width	Front wheel diameter	Rear wheel diameter	Width folded	Total width	Distance between handles	Handle heights from medium to large
Large 9.2 kg	up to 200 kg	Large 62 cm	56 cm	285 x 44 mm	235 x 44 mm	25 cm	69 cm	55 cm	77-112 cm
Medium 9.0 kg excl. shopping bag 390 grams		Medium 53 cm							

BESONDERHEITEN



Large front wheels ideal for obstacles




Including cane holder




Ergonomic grips moulded to shape



Brake lever can be operated without much effort



Ergonomic handle grips for lifting and carrying



Proven and secure "Click & Safe" locking system

ACCESSORIES



Back strap



LED-Lamps



Comfort seat



Your Rehasense Dealer

Rehasense's quality management complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TÜV NORD. Rehasense's medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR).

