



ICON 35

LIGHT & COMPACT



ICON 35

THE TRANSPORT WHEELCHAIR SERIES



from
8.2 kg



Large size
RWL35D6IB0 (Lx)
RWB35D6IB0 (Bx)

Medium size
RWL35D4IB0 (Lx)
RWB35D4IB0 (Bx)

Ice Blue



NEW: OUR TRANSPORT WHEELCHAIRS ICON 35 LX AND BX

The ICON 35 Lx which weighs only 9.0 kg is fitted with enhanced ergonomic seating, back upholstery and adjustable arm rests for increased levels of support and pressure relief if required.

The Icon 35 Bx weighs only 8.2 kg and is the lightest wheelchair in our ICON range. This wheelchair features an ergonomic seat and back angle to provide a superb seating position.


Both models, Lx and Bx, have seat widths of 16"/40cm (18" between armrests) and 18"/45cm (20" between armrests). For transportation and storage both models can be easily folded. Taking up little space with length, height and width dimensions of just 26"x29"x11".

With a maximum user weight of 125 kg both models offer a great lightweight wheelchair solution for most users. Standard features on the ICON 35 Bx and Lx include swing away leg rests with release lever, stepwise brake control with cable, armrest (height adjustable on the LX) and carrying handle for transporting when folded.

TECHNICAL DETAILS

Weights from:	Max. user weight	Front seat height	Rear seat height	Width folded	Height folded	Lenght folded	Total width	Total heighth	Total length	Seat depth	Seat widths	Seat widths between armrests	Arm rest heights
9.0 kg (Lx)	125 kg	44-49 cm	44 cm	27 cm	74 cm	65 cm	61 cm	91 cm	105 cm	40 cm	40 cm (Lx, Bx)	45 cm (Lx, Bx)	18-25 cm (Lx)
8,2 kg (Bx)											45 cm (Lx, Bx)	50 cm (Lx, Bx)	20 cm (Bx)

SPECIAL FEATURES



User-friendly, stepwise brake control with cable

Swing away leg rest with release lever

Clever carrying function

Ergonomic back upholstery (Lx, Bx). Lx has a high-comfort padding

View of the armrests Lx and Bx. The arm-rests on the Lx can also be adjusted in height

Thanks to the „Click & Safe“ the armrests are locked in place with easy adjustment



Lateral offset on back cane for ergonomic operation

14" PU rear wheel incl. step pedal

7" PU front wheel with alu fork and glass-fibre reinforced plastic bearing house



Your Rehasense Dealer

Rehasense's quality management complies with both the ISO 9001 and ISO 13485 standards, as confirmed by TÜV NORD. Rehasense's medical devices are CE-marked and meet the requirements of EU Regulation 2017/745 (MDR).

